



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
Abbott	Lifecare 5100, Lifecare 5200	3	No			EG11481
Abbott	Lifecare 5100, Lifecare 5200	5	No			EG11482
Abbott	AAEM, AAXM, ABBV, ABEM, ABTM, ABBY series, ARTM, DHS, EK-40, EK-50, EKP-1, EKR50, EO-8, EX50	3	Yes			EG11483
Abbott	AAEM, AAXM, ABBV, ABEM, ABTM, ABBY series, ARTM, DHS, EK-40, EK-50, EKP-1, EKR50, EO-8, EX50	5	Yes			EG11484
Biochem	3101, 6100, 6200, 9100	3	Yes			EG11483
Biochem	3101, 6100, 6200, 9100	5	Yes			EG11484
Burdick	Medic 3, 4, 5	3	Yes			EG11483
Burdick	Medic 3, 4, 5	5	Yes			EG11484
Colin Medical	BP308, BP408, CBM3000	3	Yes			EG11483
Colin Medical	BP308, BP408, CBM3000	5	Yes			EG11484
Criticare	504US	3	No			EG11481
Criticare	504US, 506, 507, 508, 1100, 1508	5	No			EG11482

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
Datascope	2000, 2001, 2002, 2100, 2200, 3000, 3000A, Balloon pump, M/D3, M/D3A, MD4, Passport, Passport 2, Passport EL, Passport XG, Passport XL, System 90, System 90T, System 97, System 98XT	3	No			EG11481
Datascope	2000, 2001, 2002, 2100, 2200, 3000, 3000A, Balloon pump, M/D3, M/D3A, MD4, Passport, Passport 2, Passport EL, Passport XG, Passport XL, System 90, System 90T, System 97, System 98XT	5	No			EG11482
Datex-Ohmeda	AS/3, AS/5	3	No			EG11479
Datex-Ohmeda	AS/3, AS/5	5	No			EG11480
Datex-Ohmeda	EXCEL 210-SE	3/5	Yes			EG11478
Drager	Vitalert 1000, 2000	3	Yes			EG11483
Drager	Vitalert 1000, 2000	5	Yes			EG11484
Epic System	ES100	3	No			EG11481
Epic System	ES100	5	No			EG11482
GE Medical Systems (Marquette)	Eagle and Dash monitors, Solar modules, Tram modules (component of Cath Lab), Datex-Ohmeda F/5 and MacLab	3/5	No			EG11475
GE Medical Systems	1500 Defib	3	No			EG11477

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
GE Medical Systems	1500 Defib	5	No			EG11482
GE (Corometrics)	556	3	Yes			EG11483
GE (Corometrics)	556	5	Yes			EG11484
GE (Critikon)	8100T, Dinamap (8700, 8710, 8720, 9700, 9710, 9720)	3	Yes			EG11483
GE (Critikon)	8100T, Dinamap (8700, 8710, 8720, 9700, 9710, 9720)	5	Yes			EG11484
Hill Rom	HRRM 71-1, 71-2, Porta-Fib (I, II, III), System 5 (New), System 6 (Neonatal)	3	Yes			EG11483
Hill Rom	HRRM 71-1, 71-2, Porta-Fib (I, II, III), System 5 (New), System 6 (Neonatal)	5	Yes			EG11484
Interspec	880	3	No			EG11481
Interspec	880	5	No			EG11482
Invivo Research	Millenia, Omni-Trak	3	No			EG11481
Invivo Research	Millenia, Omni-Trak	5	No			EG11482
IVY Biomedical	101, 101B, 101NR, 101R, 101T, 201, 202, 301, 302, 303, 401, 404, 405, 703, 703T, 702A	3	No			EG11481
IVY Biomedical	101, 101B, 101NR, 101R, 101T, 201, 202, 301, 302, 303, 401, 404, 405, 703, 703T, 702A	5	No			EG11482

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

Fetal Monitoring

Pulse Oximetry

Cardiac Monitoring

Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
MDE	20100, 20300, 400673, Escort (100, 101, 102, 200, 300), Prism	3	Yes			EG11483
MDE	20100, 20300, 400673, Escort (100, 101, 102, 200, 300), Prism	5	Yes			EG11484
Medtronic Physio Control	VSM (1, 3, 4, 5)	3	Yes			EG11483
Medtronic Physio Control	VSM (1, 3, 4, 5)	5	Yes			EG11484
Medtronic Physio Control	Lifepak 6, Physio 1440/R, VSM 2 (without resistors)	3	No			EG11481
Medtronic Physio Control	911ST, Lifepak (5, 6, 6S, 7, 8, 9, 9A, 9P, 10), Physio 1440/R, VSM2	5	No			EG11482
Medtronic Physio Control	911ST, Lifepak (5, 6S, 7, 8, 9, 9A, 9P, 10)	3	No			EG11477
Motorola	Apcor Radio	3	No			EG11477
Motorola	Apcor Radio	5	No			EG11482
Narco Scientific	Porta-Fib 2, 3	3	Yes			EG11483
Narco Scientific	Porta-Fib 2, 3	5	Yes			EG11484
Nellcor	OEC 6120A, N200, N250	3	Yes			EG11483
Nellcor	OEC 6120A, N200, N250	5	Yes			EG11484

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

Fetal Monitoring

Pulse Oximetry

Cardiac Monitoring

Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



















Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
Nihon Kohden	BSM Series, Cardiolife 5000 series, 7000 series, 8000 series defibrillator, Lifescope 6, 8, 11, 12, OEC 6102A (Serial#>20156), OEC 6105A (Serial#>20226), OEC 7102A, OEC8108A, TEC 5200A (Serial#>20421), TEC 7100A, TEC 7200A, TEC 7300A	3	Yes			EG11483
Nihon Kohden	BSM Series, Cardiolife 5000 series, 7000 series, 8000 series defibrillator, Lifescope 6, 8, 11, 12, OEC 6102A (Serial#>20156), OEC 6105A (Serial#>20226), OEC 7102A, OEC8108A, TEC 5200A (Serial#>20421), TEC 7100A, TEC 7200A, TEC 7300A	5	Yes			EG11484
Nihon Kohden	Life Scope 9	3	No			EG11481
Nihon Kohden	Life Scope 9	5	No			EG11482
Novametrics	Life Defense I	3	Yes			EG11483
Novametrics	Life Defense I	5	Yes			EG11484
Novametrics	Neonatal 902, 903, 910	3	No			EG11481
Novametrics	Neonatal 902, 903, 910	5	No			EG11482

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

— Fetal Monitoring

— Pulse Oximetry

— Cardiac Monitoring

— Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
Philips	78352C, 78354C, 78834C Neonatal Bedside, Codemaster, M1001A/B, M1002A/B, M1106A, M1175A, M1176A, M1280A, M1722A/B, M1723A/B, M1500A, M1520A, M1540C, M1560C, M1723A/B, M2475B, M3000A, M3001A, M4735A, Merlin System, OmniCare	3/5	Yes			EG11478
Philips (ATL)	800, 800 Plus, Apogee, CX, CX200, CX400, HDI 9, HDI 3000, HDI 5000, Interspec Cardiology XL, Mark (IV, VIII), Ultramark (4, 4PLUS, 5, 7)	3	Yes			EG11483
Philips (ATL)	800, 800 Plus, Apogee, CX, CX200, CX400, HDI 9, HDI 3000, HDI 5000, Interspec Cardiology XL, Mark (IV, VIII), Ultramark (4, 4PLUS, 5, 7)	5	Yes			EG11484
Pioneer	26125, 33167, 5260, 5330, Pulsar 4	3	Yes			EG11483
Pioneer	26125, 33167, 5260, 5330, Pulsar 4	5	Yes			EG11484
PPG (Datamedix)	Guardian 1001, N1001	3	Yes			EG11483
PPG (Datamedix)	Guardian 1001, N1002	5	Yes			EG11484
Siemens Medical Systems	Sirecust 300DU	3	No			EG11481
Siemens Medical Systems	Sirecust 300DU, 1200 series, 400 series, 500 series, 700 series, 900 series	5	No			EG11482
Spacelabs Medical	900 series with respiration	3	No			EG11481

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
Spacelabs Medical	900 series with respiration, 90601A	5	No			EG11482
Spacelabs Medical	500 Series, 600 Series, 900 Series, 90300, 90308, 90308-11, 90303A/B, 90305, 90418, 90431, 90470-01, 90470-02, 90600, 90602/A, 90603/A, 90604A, 90621A, 90622A, 90623A, 90624A, 90662A, 90700, 90701, 90721, 90724, PC Express, PC, PCMS, PC II	3	Yes			EG11483
Spacelabs Medical	500 Series, 600 Series, 900 Series, 90300, 90308, 90308-11, 90303A/B, 90305, 90418, 90431, 90470-01, 90470-02, 90600, 90602/A, 90603/A, 90604A, 90621A, 90622A, 90623A, 90624A, 90662A, 90700, 90701, 90721, 90724, PC Express, PC, PCMS, PC II	5	Yes			EG11484
Spacelabs Medical	90601A	3	No			EG11477
Spacelabs Medical	90385 with Command Module 90496	3/5	Yes			EG11476
Spacelabs (Vitatek)	4014	3	No			EG11481
Spacelabs (Vitatek)	4014	5	No			EG11482
Spacelabs (Vitatek)	208 (with Lead Select), 4045 Defib, 408(with Lead Select), 412 (with Lead Select), 413 (with Lead Select), 414 (with Lead Select), 511, 512, 514, 521, 522, 524, 526	3	Yes			EG11483

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

— Fetal Monitoring

— Pulse Oximetry

— Cardiac Monitoring

— Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337



Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

All Patient Cables and Leadwires have received FDA 510(k) clearance

Monitor OEM	Monitors Supported	Number of Leads	1K Ω Resistor	Monitor Connector	Trunk Cable Connector	AIV Part Number
Spacelabs (Vitatek)	208 (with Lead Select), 4045 Defib, 408(with Lead Select), 412 (with Lead Select), 413 (with Lead Select), 414 (with Lead Select), 511, 512, 514, 521, 522, 524, 526	5	Yes			EG11484
Temtech	Life Guard 9, M401, M501, M503	3	Yes			EG11483
Temtech	Life Guard 9, M401, M501, M503	5	Yes			EG11484
Welch Allyn (MRL)	360SL, 450SL, 450 PreAmp, 521A, PortaPac 80, 90, 1000	3	Yes			EG11483
Welch Allyn (MRL)	360SL, 450SL, 450 PreAmp, 521A, PortaPac 80, 90, 1000	5	Yes			EG11484
Welch Allyn (Protocol)	Propaq 100 Series, 200 Series, Propaq Encore, Ultra Smart	3	Yes			EG11483
Welch Allyn (Protocol)	Propaq 100 Series, 200 Series, Propaq Encore, Ultra Smart	5	Yes			EG11484
Zoll	D900, M Series, PD900 PD1200, PD1400	3	Yes			EG11483
Zoll	D900, M Series, PD900 PD1200, PD1400	5	Yes			EG11484

AIV Patient Cables must be used in combination with AIV Leadwire Sets for the system to function as intended.

(continued)

Call us today 800.990.2911 or visit us at www.aiv-inc.com

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

Specialized Repair Services
Accessories
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337




Cardiac Monitoring

ECG Cables and Leadwires

Cardiac Monitoring

- **Reliable signal** - Fully shielded patient cables and leadwires include braided metal shielding to reduce RFI and EMI
- **Dependable performance** - Engineered to be functionally equivalent to OEM models
- **Convenient ordering** - All items are shipped from stock for fast delivery
- **Protection for your investment** - Cables are covered by a 1-year warranty and leadwires by a 30-day warranty

All Patient Cables and Leadwires have received FDA 510(k) clearance

Type of Lead	Number of Leads	AIV Patient Cable Supported	130M Ω Resistor	Leadwire Connector	Electrode Connector	AIV Part Number
Snap	3	EG11475	Yes			EG11697
Snap	3	EG11476, EG11477, EG11478, EG11479, EG11481, EG11483	No			EG11487
Snap	5	EG11475, EG11476, EG11478, EG11480, EG11482, EG11484	No			EG11485
Grabber	3	EG11475	Yes			EG11706
Grabber	3	EG11476, EG11477, EG11478, EG11479, EG11481, EG11483	No			EG11488
Grabber	5	EG11475, EG11476, EG11478, EG11480, EG11482, EG11484	No			EG11486

Specialized Repair Services
Accessories
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337

**Call us today 800.990.2911 or
visit us at www.aiv-inc.com**

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.