

IsoMATE™



For Use with Medical Equipment in Patient Care Areas

- Fully contained isolation transformer offers full line isolation
- Less than 100 μ A leakage current for all connected equipment
- Dual circuit breakers provide overload protection
- Regulatory Approvals: UL listed to 60601-1, CSA C22.2 No.601
- Mounts to an IV pole with a diameter of 0.75" to 2.50"
- Design features include built-in handle and locking hex nut. Angled surface allows for connection of 3 wall mount AC adapters and in-line plugs simultaneously.
- 2-year warranty



Call us at 800.990.2911 or visit www.aiv-inc.com to learn more.

AIV is an ISO 13485 registered facility.