













# Pulse Oximetry Adapter Cables

AIV offers a large line of adapter cables for pulse oximeters and pulse oximetry modules, including the latest generation of technology from Nellcor® and Philips®. These cables are functionally equivalent to the OEM and contain upgraded shielding to reduce EMI and enhance signal quality.

Items are shipped from stock • 90-day warranty

## Adapter Cables

Monitor OEM	Devices Supported	SpO2 Technology	Monitor Connector	Sensor Connector	OEM Part Number	AIV Part Number
Datascope®	Passport®, Sat 3, Sat 4	Datascope®			0012-00-0516-02	CBDAT CE
GE® Medical Systems (Marquette®)	GE Tram Module T451N and T851N with OxiSmart technology, TRAM 2250, TRAM 650, TRAM 850, DASH, Eagle, Solar, Mac-Lab Cardiac Catheterization System	Nellcor®			E90004GE 2000-6644-0001	CB10640 CE
Nellcor®	NPB-290, N-295, N-395	Nellcor®			MC10	CB10587 CE
Nellcor®	Nellcor N-20, NPB-40, NPB-75, NPB-190 Welch Allyn® Protocol Propaq®, Most BCI monitors	Nellcor®			DEC8	CB10588 CE
Nellcor®	Nellcor N-595, GE Dinamap® Pro Series with OxiMax® Modules	Nellcor®			DOC10	CB10639

(continued on back)

Specialized Repair Services  
Accessories  
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.  
7485 Shipley Avenue  
Harmans, MD 21077 USA  
p 410.787.1300  
f 410.787.1337

**Call us today 800.990.2911 or  
visit us at [www.aiv-inc.com](http://www.aiv-inc.com)**

\*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

















# Pulse Oximetry Adapter Cables

AIV offers a large line of adapter cables for pulse oximeters and pulse oximetry modules, including the latest generation of technology from Nellcor® and Philips®. These cables are functionally equivalent to the OEM and contain upgraded shielding to reduce EMI and enhance signal quality.

**Items are shipped from stock · 90-day warranty**

## Adapter Cables

(continued from front)

Monitor OEM	Devices Supported	SpO2 Technology	Monitor Connector	Sensor Connector	OEM Part Number	AIV Part Number
Nellcor®	NPB-290, N-295, N-395, N-3000	Nellcor®			SCP10	CB10641 CE
Philips®	M1175A (with module M1020A), M1176A (with module M1020A), 78352C, 78354C, M1025B, M1205A, M1722A, M1722B, M1723B, M2475B	Nellcor®			M-1900-B	CBHPN1
Philips®	M1175A (with module M1020A), M1176A (with module M1020A), 78352A, 78352C, 78354A, 78354C, 78833A, M78833B, M1722A, M1722B, M1723A, M1723B, M2475B	Philips®			M-1940-A	CBHPHP1 CE
Philips®	M3000A-Viridia® Multi-Measurement Monitor, M3500B-Heartstream® XLT Defib/Monitor, M4735A-Philips® Defib/Monitor	Philips®			M1914A	CB10586 CE
Philips®	M2601A-HP Viridia® Patient Monitor, M3000A-Viridia® Multi-measurement Monitor, M3500B-Heartstream® XLT Defib/Monitor, M4735-A Philips® XL Defib/Monitor	Nellcor®			M1943A	CB10590 CE
Spacelabs®	Monitors with SpO2 modules: 90351-D/G, 90465/66/67, 90489, 90651A-08	Nellcor®			700-0002-00	CBS P1 CE
Spacelabs®	Ultraview® Series 1050, 1600, 1700 All monitors using the Ultraview® Command Module	Nellcor®			700-0030-00	CBS P2 CE

Specialized Repair Services  
Accessories  
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

**Call us today 800.990.2911 or  
visit us at [www.aiv-inc.com](http://www.aiv-inc.com)**

\*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.

AIV, Inc.  
7485 Shipley Avenue  
Harmans, MD 21077 USA  
p 410.787.1300  
f 410.787.1337