



We Extend Your Life!

AIV Parts & Repairs for Baxter Flo-Gard 6201

**AIV Parts & Repairs For
Baxter Flo-Gard 6201**

Everything you need to keep your Baxter Flo-Gard 6201 pumps alive.

- 6 month warranty on repairs
- 7 to 10 day turnaround on repairs
- 120 day warranty on new replacement parts
- 60 day warranty on recertified parts
- Service agreements and rentals available - call for details
- All device repairs completed by trained technicians in our ISO 13485:2016 certified facility



AIV New Replacement Parts for the Baxter Flo-Gard 6201		
Description	AIV Part Number	OEM Reference Number
Door Assembly Right (P2)	BX10468	F049120021
Door Latch Assembly Right (P2)	BX10675	CHNDP100.F
Front Case 6201	BX10610	GCABA104A
Keypad 6201	HPNL1025A	HPNL1025.A
Label Set Complete 6201	BX11073	072601081
Cover Pump Head 6201	BX13029	GCOVA101A
Power Cord 8 foot	PCM8	QACCD761.D
Rear Case 6201	BX10938	GCABB104C
Tube Misloading Sensor L-P2	10666	QSWZ1011A
Tube Misloading Sensor R-P2	10667	QSWZ1012A
Thumb Shaft Assembly	BX12828	N/A
AIV Certified Parts for the Baxter Flo-Gard 6201		
Description	AIV Part Number	OEM Reference Number
Assembly Air Sensor Pump 2 6201	S048993RC	F049120023
Assembly Pole Clamp Flo-Gard	BX12685RC	F04912008
Backup Audio Alarm Flo-Gard	CPWBN104CRC	CPWBN104C
Battery Fuse 2A 125V Flo-Gard	QFSB4451ARC	QFSB4451A
Battery Harness 6201	QCNW1094ARC	QCNW1094A
Buzzer Harness 6201	QCNW1093ARC	QCNW1093A
Cable Flat Ribbon Flo-Gard	QCNW1117ARC	QCNW1117A
Cover Battery Compartment 6201	GFTAB101ARC	GFTAB101A
Cover Power Cord Flo-Gard	GCOVH101BRC	GCOVH101B
CPU Board 6201	CPWBX104ARC	CPWBX104A
Display Board 6201	CPWBN104ARC	CPWBN104A
Knob Panel Lock Switch Flo-Gard	JKNBZ101CRC	JKNBZ101C
Knob Volume Control Flo-Gard	JKNBZ100CRC	JKNBZ100C
Occlusion sensor Flo-Gard	OCCSENI1AIVRC	N/A
PCB assembly 6201	CPWBN105ARC	CPWBN105A
PCB Sensor Board 6201	CPWBX104BRC	CPWBX104B
Power Cord Strap Flo-Gard	UBNDA100ARC	UBNDA100A
Power Receptacle Flo-Gard	QSOCA101ARC	QSOCA101A
Power supply PCB 6201	RDENC100DRC	RDENC100D
Safety Slide Clamp assembly Flo-Gard	CCRA1006ARC	CCRA1006A
Seal Pump Head Flogard	PGUMM103KRC	PGUMM103K
Seal Rubber housing 6201	PGUMM103JRC	PGUMM103J
Terminal PCB Flo-Gard	CPWBN104BRC	CPWBN104B
Transformer Power 6201	RTRNP100ERC	RTRNP100E
Vol / Lock Ground Plate 6201	LANGQ103BRC	LANGQ103B
Other parts available - Call for details		
AIV Repair Services for the Baxter Flo-Gard 6201		
Description	AIV Part Number	OEM Reference Number
6201 Repair - fixed price labor (parts are additional)	RP6201	N/A
6201 Flat Rate Repair - price includes labor and any necessary parts	FLRT6201	N/A
CPU Printed Circuit Board Repair	PBCPWXB104A	CPWBX104A
Sensor Printed Circuit Board Repair	PBCPWXB104B	CPWBX104B



Call us today
877.266.6897
or visit us at
www.aiv-inc.com

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 877.266.6897
f 410.787.1337
aivsales@aiv-inc.com



The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.