



America

CERTIFICATE

No. U8 13 05 84362 001

Holder of Certificate: American I.V. Products, Inc.
dba AIV, Inc.
7485 Shipley Avenue
Harmans MD 21077
USA

Production Facility(ies): 84525

Certification Mark:



Product: General Medical Devices
Relocateable Power Tap

Model(s): Powermate 60601 PS12944

Parameters:

Input voltage:	120VAC
Input frequency:	50/60Hz
Input Current:	15A
Type of protection against electric shock:	Class I
Degree of protection against electric shock:	N/A – No applied parts
Degree of protection against ingress of liquids:	IP21
Protection against flammable anesthetics:	Not protected
Mode of operation:	Continuous

Tested according to: UL 60601-1:2003
CAN/CSA-C22.2 No. 60601-1-08
IEC 60601-1/A2:1995
IEC 60601-1:2005

The product was voluntarily tested according to the relevant safety requirements noted above. It can be marked with the certification mark above. The mark must not be altered in anyway. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC Guide 67. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited certification body.

Test report no.: DM1211737

Date: 2013-05-24

TÜV SÜD AMERICA INC • 10 Centennial Drive • Peabody, MA 01960 USA • www.TUVamerica.com



TUV®